



CHIRAG MEDITECH PVT. LTD.

Chirag Meditech Pvt Ltd is the leading manufacturer of Trans Urethral Resection of Prostrate (TURP) electrodes, used for surgical treatment of enlarged prostrate / BPH, for Domestic and Export markets.

Electrodes are of different types for various purposes:

Cutting Electrode : Most commonly used electrode, removes the prostate tissue in the form of small chips.


Colling's Electrode: is used to make straight incisions in very hard prostate tissue.

Roller Electrode: used to vaporize the prostate tissue using high power.








Ball Electrode: is used at the end of surgery, some surgeons use this electrode to coagulate and smoothen the freshly resected surface of the prostate gland.

Cold Knife : The strictures can be incised with a cold knife.

Chirag Meditech Pvt Ltd manufactures different design / make of Loops compatible to resectoscopes of Karl Storz, Richard Wolf, Olympus, etc.



Chirag Loop Electrodes are manufactured with high degree of precision with R&D as main focus and are made to pass stringent quality tests.

Design / Make							
	Cutting	Collings Knife	Roller	Small Ball - 3mm	Big Ball - 5mm	Cold Knife	Cutting 27F
Storz Double Stem 30° - 24F	CDC24F	CDZ24F	CDR24F	CDS24F	CDB24F	CDK24F	CDC27F
Storz Single Stem 30° - 24F	CSC24F	CSZ24F	CSR24F	CSS24F	CSB24F	CSK24F	CSC27F
Richard Wolf Single Stem 15° - 24F	CWC24F	CWZ24F	CWR24F	CWS24F	CWB24F	CWK24F	CWC27F
Olympus Single Stem 12° - 24F	COC24F	COZ24F	COR24F	COS24F	COB24F	COK24F	-
Olympus Single Stem 30° - 24F	CLC24F	CLZ24F	CLR24F	CLS24F	CLB24F	CLK24F	-

Quality : Chirag Turp Loop electrodes are High Voltage Tested for insulation and practically tested with various settings using both Electro surgical Unit and Diathermy machine and working element unlike other competitive manufacturers.

We can also develop TURP Loop electrodes to any designs / make of resectoscopes.



ISO 9001-2008, ISO 13485-2003 & CE

Regd office: #2234, 9th main, 'E' Block, Rajajinagar 2nd Stage, Bangalore - 560 010, INDIA.

Factory: No.28/1, 1st floor, 4th Main, Industrial Town Rajajinagar, Bangalore - 560044, INDIA

Mobile No. : +91 98451 26320, +91 94498 26320, Fax No: +91 80 23157783

Email: chiragmeditech@yahoo.com, info@chiragmeditech.com, chiragmeditech@hotmail.com

MFRS OF TRANS URETHRAL RESECTION OF PROSTATE ELECTRODES & UROLOGY PRODUCTS